

## AUTOCLAVES



Portable / Vertical

MOC : Full SS

Operating Pressure: 15 – 18 psi

Operating Temp.: 3 Heat Control Switch

SIZE	Weight	Vol
12" Ø X 12" Ht	11 kgs	15 ltr
12" Ø X 20" Ht	37 kgs	35 ltr
14" Ø X 22" Ht	45 kgs	52 ltr
18" Ø X 24" Ht	60 kgs	95 ltr
22" Ø X 30" Ht	80 kgs	178 ltr

Accessories: Baskets / Racks

Optional: Pressure Control Switch, Low Water Level Cut Out, Digital temp Controller, Timer, Vacuum Breaker Valve