

PORTABLE AUTOCLAVES

APPLICATION :

- Clinics / Nursing Homes / Hospitals / Medical Consultants
- Laboratories/ Medical Waste Treatment Facilities / Animal Care Facilities
- Pharmaceutical / Biotechnology & Allied Industries

CONSTRUCTION / FEATURES : Single walled All SS Construction with all joints argon arc welded and grounded to crevice free finish, The Top cover & support Flange made of thick SS Plate die pressed to dome shaped and machined. Electrically operated Operating Pressure / temp. 15Psi/121°C

STANDARD FITTINGS :

- Pressure Gauge, Safety Valves, Steam Release Valve, Neoprene/Silicon Gasket

OPTIONAL : Digital Temperature Controller cum indicator / Solenoid valve.

